90295632

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
Silk · expert Pro

## How does it work?

Silk-expert Pro is designed to help break the cycle of hair growth. Light energy is transferred through the skin's surface and is absorbed by melanin present in the hair shaft. The absorbed light energy is converted to heat energy (below the surface of the skin), which disables the hair follicle and puts it to sleep.

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## Skin tone and hair color check

Refer to the chart to check whether your skin tone and hair color are suitable for IPL treatment.

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## Skin Pro 2.0 (SensoAdapt™)

Continuous and automatic adjustment of light intensity – 80x per second.

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- < **GLIDING METHOD** for large areas:  
Press treatment button and hold while gliding along the leg.

## How to treat


- STAMPING METHOD** for small areas:  
Press and release the treatment button, lift off the treatment head and place it on the next adjacent spot. Repeat spot by spot.

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
## < Regimen

1. Start-up phase: once per week for the first 4 to 12 weeks
2. Maintenance: continue to your needs (e.g. every 1–2 months)

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## < Prepare

Remove visible hairs, clean and dry the area you want to treat. Perform a patch test and wait 24 hours before the first full treatment.

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
## Start

Connect the device to an electric outlet. Press any button to activate the device.

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## < How to flash

1. White LEDs: valid skin tone is detected
2. No light: no skin contact
3. Red LED: invalid skin tone / no full contact

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## SensoAdapt™ SKIN SENSOR

Continuously & automatically  
reads your skin tone and  
adapts the energy for safe and  
effective use

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### TREATMENT BUTTON

Hold for fast gliding  
treatment or press & release  
for precise treatment

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
### SENSITIVE & EXTRA SENSITIVE MODE

For first-time users and sensitive  
areas

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More details on:

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Braun.com

**PRECAUTIONS:** DO NOT USE the device on any area where you may want hair to grow back. Hair removal by lasers or intense pulse light sources can cause increased hair growth in some individuals.

**WARNINGS:** DO NOT use if you have a known sensitivity to sunlight (photosensitivity) or are taking a medication that causes photosensitivity. Always check the instructions that come with your medical products to see if photosensitivity is listed as a side effect. Using the device on photosensitive skin could result in skin damage including blisters and swelling.

**RISKS OF USING THE SILK-EXPERT PRO:** The side effects reported in the clinical trial were: Itchiness, Slight Redness, Mild Stinging Sensation, In-growing Hair, Follicle Redness. In all cases, the treated skin returned to normal within 7 days.

**CONTRAINDICATIONS:** DO NOT use the device if you have very dark skin. DO NOT use the device if your skin has been artificially or naturally tanned or is burned from sun overexposure. DO NOT expose treated areas to the sun. Wait at least 7 days after treatment before exposing treated skin to direct sunlight. DO NOT use the device directly on the nipples, genitals or around the anus. DO NOT use the device if you have a history of skin cancer or pre-cancerous lesions. DO NOT use if you are pregnant, lactating. DO NOT use on dark brown or black spots such as birthmarks, moles or warts. DO NOT use if you have any skin condition in the treatment area, including psoriasis, vitiligo, eczema, acne, herpes simplex or active infections or wounds. DO NOT use if you are under the age of 18. DO NOT use on tattoo or permanent make-up in the treatment area. DO NOT use on beard/facial hair for men.

**WARNINGS:** The Silk-expert Pro is an electrical device that emits intense pulses of light. DO NOT look directly into the treatment window.

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